



## **Sensores de temperatura de pastillas de freno de termopar IEC tipo K - 4,0 mm de diametro**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Estos sensores de termopar se utilizan para controlar la temperatura de las pastillas de freno hasta +400 C. Estn diseados para incrustarse en las almohadillas.

Cada sensor tiene 1 m de cable aislado de fibra de vidrio, con trenzado de acero inoxidable. Disponible en termopar tipo K o J y con una seleccin de dimetros de punta. Cada punta est sellada con una pasta cermica de alta temperatura (Fortafix) que puede soportar temperaturas de hasta +400 C.

### Caractersticas tcnicas

## Specifications

|                       |                                                       |
|-----------------------|-------------------------------------------------------|
| Product Code          | XE-3894-001                                           |
| General Description   | Designed for monitoring the temperature of brake pads |
| Sensor Type           | Grounded Junction                                     |
| Thermocouple Type     | K                                                     |
| Probe Diameter        | 4.0mm                                                 |
| Probe Length          | Tip Length: 4mm                                       |
| Cable Length          | 1M                                                    |
| Cable Type            | Fibre glass with stainless steel Overbraid 7/0.2mm    |
| Termination Type      | Bare Tails                                            |
| Max. Temperature      | 400°C                                                 |
| Standards Met         | IEC                                                   |
| Important Information | Potting sealant: Ceramic paste (Fortafix)             |