


Conector de termopar de montaje en miniatura con soporte de acero inoxidable FMTC-CU-SSPF Tipo Cu



Miniature in-line Socket with stainless steel mounting bracket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar de montaje en miniatura con soporte de acero inoxidable tipo FMTC-CU-SSPF

Soporte de acero inoxidable simplemente montado en el conector. Terminales totalmente cerrados con sellado de puerto de cable

Especificaciones**Specifications**

Product Code	XE-1378-001
General Description	Fully enclosed terminals with cable port seal
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	Socket with SS Bracket
Connector Size	Miniature
Colour	White
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS