



**LABFACILITY**  
TEMPERATURE & PROCESS TECHNOLOGY

## **Termopar con aislamiento mineral de respuesta rápida de 0,25 mm de diámetro con enchufe en miniatura ANSI - Tipo K**



IMAGE  
COMING  
SOON

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockists of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Utilice nuestros [formularios de sensores personalizados](#) para enviar sus requisitos para varias longitudes de MI, ya que todas terminan con un conector de 220 C para una conexión rápida y sencilla.

### **Necesita calibrar su sensor de temperatura?**

Trazable a los estándares de calibración ISO17025 (UKAS), Labfacility ofrece la opción de una calibración trazable de 3 o 5 puntos utilizando nuestras instalaciones de calibración internas.

**Tenga en cuenta que si agrega una opción de calibración a su pedido, se comunicará con usted con respecto a la información de entrega.**

### **Características técnicas**

## **Specifications**

<b>General Description</b>	Highly flexible, sheath can be bent/formed to suit many applications and processes
<b>Sensor Type</b>	Type K ANSI MB STYLE probes 0.25mm
<b>Thermocouple Junction</b>	Single element, junction insulated or grounded
<b>Probe Length</b>	100mm / 150mm / 250mm
<b>Sheath Material</b>	Stainless steel 316 sheath
<b>Termination Type</b>	Fitted Miniature Plug Type ANSI
<b>Range</b>	-40 to +260°C
<b>Important Information</b>	A 1%+ Length Tolerance applies