

Sonda termistor calibrable Lascar EL-P-TPX+



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sonda termistor calibrable de alta precisin, de 3 metros de longitud con punta de acero inoxidable, para su uso con EL-WiFi-TPX

[+ \(XE-4076-001\)](#)

Especificaciones

Specifications

Product Code	XE-4081-001
Model	EL-P-TPX+
Type	Calibratable thermistor probe
Accuracy	+/-0.2°C (-15 to +80°C) / +/-0.6°C (-40 to 125°C)
Cable Type	3m (118.1")
Measuring range	-40 to +125°C / -40 to + 257°F