

Sonda termistor de gliclic Lascar EL-P-VACX de alta precisin, sin calibrar



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sonda termistor de 3 metros de longitud y de alta precisión, sin calibrar, en un frasco sellado lleno de glicol para el monitoreo de vacunas. Para uso con el registrador de datos EL-WiFi-VACX

[\(XE-4058-001\)](#)

Especificaciones

Specifications

Product Code	XE-4082-001
Model	EL-P-VACX
Type	High accuracy Glycol thermistor probe, uncalibrated
Accuracy	+/-0.1°C / +/-0.2°F
Cable Length	3 metres
Measuring range	-40 to +60°C / -40 to + 140°F