




Conector de termopar para montaje de PCB IM-N-PCB Tipo N IEC



Miniature PCB Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar para montaje de PCB IM-N-PCB Tipo N IEC

Diseado especialmente para PCB, pines slidos para montaje directo en placa y provisin para sensor CJC

Especificaciones**Specifications**

Product Code	XE-1259-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	N IEC
Connector Type	PCB socket
Connector Size	Miniature
Colour	Pink
Max. Temperature	220°C
Pins	Brass
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS