


Conector de termopar para montaje de PCB Socket FMTC-CU-PCB Type Cu



Miniature PCB Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar para montaje de PCB Socket FMTC-CU-PCB Tipo CU

Diseado especialmente para PCB, pines slidos para montaje directo en placa y provisin para sensor CJC

Especificaciones**Specifications**

Product Code	XE-1379-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	PCB socket
Connector Size	Miniature
Colour	White
Max. Temperature	220°C
Pins	Brass
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector