



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar de respone rpida de 0,5 mm de dimetro con olla roscada y colas IEC Tipo K - Unin conectada a tierra - 250 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Tipo K con olla roscada y colas de acero inoxidable 310

Para las versiones de 0,25 mm, 0,75 mm y 0,81 mm de diámetro y para las versiones de funda Inconel Alloy 600, utilice nuestros

[formularios de sensores personalizados](#) para enviar sus requisitos o envíe un correo electrónico a sales@labfacility.com / llame al 871280

Necesita calibrar su sensor de temperatura?

Trazable a los estándares de calibración ISO17025 (UKAS), Labfacility ofrece la opción de una calibración trazable de 3 o 5 puntos utilizando nuestra instalación de calibración interna, puede seleccionar cualquier punto de temperatura entre -10 °C y +200 °C.

Tenga en cuenta que si agrega una opción de calibración a su pedido, nos pondremos en contacto con usted con respecto a la información de entrega.

Características técnicas

Specifications

Product Code	XE-3292-001
General Description	Type K with threaded pot & tails in 310 stainless steel
Sensor Type	Mineral Insulated
Thermocouple Type	K
Element/hot junction	Single element, junction grounded from sheath
Probe Diameter	0.5mm
Probe Length	250mm
Sheath Material	310 stainless steel
Pot Type	M8 x 1.0mm fine pitch threaded pot seal (200°C)
Cable Length	100mm
Cable Type	PFA insulated twin twisted IEC
Core / Strands	7/0.2mm



Termination Type	Bare Tails
Max. Temperature	+650°C
Min. Temperature	-40°C
Tolerance	Class 1
Standards Met	In accordance with IEC 584