


**Sello liso de olla de 0,81 mm de diametro con 1 metro de plomo aislado con PFA, par plano ANSI 7/0,2 mm con colas Unin aislada - 150 mm de longitud de la sonda**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sello de trampa liso con 1 metro de plomo aislado con PFA, par plano de 7/0,2 mm con cola.

Si necesita cantidades mayores, por favor envíe un correo electrónico a [sales@labfacility.com](mailto:sales@labfacility.com) o llame a nuestra línea de ventas al 01243 871280.

### **Necesitas calibrar tu sensor de temperatura?**

Rastreable a estándares de calibración ISO17025 (UKAS), Labfacility ofrece la opción de calibración trazable en 3 o 5 puntos utilizando nuestra instalación de calibración interna.

Ten en cuenta que si añades una opción de calibración a tu pedido, te contactaremos para información de entrega.

### **Características técnicas**

## **Specifications**

|                             |                                                                                                   |
|-----------------------------|---------------------------------------------------------------------------------------------------|
| <b>Product Code</b>         | XE-0259-001                                                                                       |
| <b>General Description</b>  | Plain pot seal (200°C) with 1 Metre of PFA insulated lead, 7/0.2mm flat pair with tails           |
| <b>Sensor Type</b>          | Mineral Insulated                                                                                 |
| <b>Thermocouple Type</b>    | K                                                                                                 |
| <b>Element/hot junction</b> | Single element, junction insulated from sheath                                                    |
| <b>Probe Diameter</b>       | 0.81mm                                                                                            |
| <b>Probe Length</b>         | 150mm                                                                                             |
| <b>Sheath Material</b>      | K - Stainless steel 316                                                                           |
| <b>Pot Type</b>             | Plain (+200°C)                                                                                    |
| <b>Cable Length</b>         | 1M                                                                                                |
| <b>Cable Type</b>           | PFA Teflon® insulated flat pair cable and tails (colour coded to ANSI MC96.1)<br>(Max Temp 250°C) |

|                  |                                |
|------------------|--------------------------------|
| Core / Strands   | 7/0.2mm                        |
| Termination Type | Bare Tails                     |
| Max. Temperature | K +650°C                       |
| Min. Temperature | K -40°C                        |
| Tolerance        | Class 1                        |
| Standards Met    | In accordance with ANSI MC96.1 |