




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar magnético de potencia de tracción (herradura) tipo K IEC de 11,8 kg, cable aislado con PFA con trenza de acero inoxidable que termina en enchufe estándar - 1,5 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-0077-001
General Description	11.8kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	45 x 30 x 30mm (L x H x W) - Max pull: 11.8 KG
Cable Length	1.5M
Cable Type	PFA Insulated Cable with Stainless Steel Over-Braid
Core / Strands	7/0.2mm
Termination Type	Standard Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC