



## Conector estndar para termopar FSTC-CU-M Tipo Cu



**Standard in-line Plug**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector estndar para termopar FSTC-CU-M Tipo CU**

Los contactos estn hechos de aleaciones verdaderas de termopar. Adecuado para aplicaciones industriales

**Especificaciones****Specifications**

<b>Product Code</b>	XE-1670-001
<b>General Description</b>	Contacts are polarised to prevent incorrect connection
<b>Thermocouple Type</b>	IEC, ANSI, JIS Type CU
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Standard
<b>Colour</b>	White
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Solid Round Pins
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper
<b>Standards Met</b>	BS EN 50212:1997. RoHS