


**Termopar magnético de botón de tracción IEC tipo K  
de 1,9 kg, cable aislado con PFA con trenza de  
acero inoxidable que termina en enchufe estándar  
- 1M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

|                            |                                                  |
|----------------------------|--------------------------------------------------|
| <b>Product Code</b>        | XE-0081-001                                      |
| <b>General Description</b> | 1.9kg Pull Button Magnet Thermocouple            |
| <b>Sensor Type</b>         | Magnet Thermocouple                              |
| <b>Thermocouple Type</b>   | K                                                |
| <b>Magnet</b>              | 19.1mm Diameter / 12.7mm height. Max Pull: 1.9kg |
| <b>Cable Length</b>        | 1M                                               |
| <b>Cable Type</b>          | PFA Flat Pair with stainless steel overbraid     |
| <b>Core / Strands</b>      | 7/0.2mm                                          |
| <b>Termination Type</b>    | Standard Plug                                    |
| <b>Max. Temperature</b>    | 250°C                                            |
| <b>Min. Temperature</b>    | -50°C                                            |
| <b>Standards Met</b>       | IEC                                              |