




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar magnético de botón de tracción IEC tipo K de 1,9 kg, cable aislado con PFA con trenza de acero inoxidable que termina en enchufe en miniatura - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caracterstias tcnicas**Specifications**

Product Code	XE-0080-001
General Description	1.9kg Pull Button Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	19.1mm Diameter / 12.7mm height. Max Pull: 1.9kg
Cable Length	2M
Cable Type	PFA Flat Pair with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC