


**Termopar magnético de potencia de tracción  
(herradura) tipo K IEC de 4,5 kg, cable aislado de  
PFA con trenza de acero inoxidable que termina  
en enchufe en miniatura - 1M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

|                            |                                                   |
|----------------------------|---------------------------------------------------|
| <b>Product Code</b>        | XE-0083-001                                       |
| <b>General Description</b> | 4.5kg Pull Power (Horseshoe) Magnet Thermocouple  |
| <b>Sensor Type</b>         | Magnet Thermocouple                               |
| <b>Thermocouple Type</b>   | K                                                 |
| <b>Magnet</b>              | 30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg       |
| <b>Cable Length</b>        | 1M                                                |
| <b>Cable Type</b>          | PFA insulated lead with stainless steel overbraid |
| <b>Core / Strands</b>      | 7/0.2mm                                           |
| <b>Termination Type</b>    | Miniature Plug                                    |
| <b>Max. Temperature</b>    | 250°C                                             |
| <b>Min. Temperature</b>    | -50°C                                             |
| <b>Standards Met</b>       | IEC                                               |