


Termopar magnético de potencia de tracción (herradura) tipo K IEC de 4,5 kg, cable aislado con PFA con trenza de acero inoxidable que termina en enchufe estándar - 1M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

Product Code	XE-0085-001
General Description	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
Cable Length	1M
Cable Type	PFA insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Standard Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC