


Conector estndar de termopar Enchufe en Inea IS-K-F Tipo K IEC



Standard in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector estndar de termopar Enchufe en Inea IS-K-F Tipo K IEC

Acepta pasadores redondos estndar, adecuados para aplicaciones industriales

Especificaciones**Specifications**

| | |
|----------------------------|---|
| Product Code | XE-1402-001 |
| General Description | Accepts standard round plug pins, contacts are polarised to prevent incorrect connection. |
| Thermocouple Type | K IEC |
| Connector Type | Socket |
| Connector Size | Standard |
| Colour | Green |
| Max. Temperature | 220°C |
| +Positive Contact | Nickel Chromium |
| -Negative Contact | Nickel Alloy |
| Standards Met | BS EN 50212:1997. RoHS |
| Series | Standard Connector |