


Conector estndar de termopar Enchufe en Inea IS-R/S-F Tipo R/S IEC



Standard in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector estndar de termopar Toma en Inea IS-N-F Tipo N IEC

Acepta pasadores redondos estndar, adecuados para aplicaciones industriales

Especificaciones**Specifications**

Product Code	XE-1592-001
General Description	Accepts standard round plug pins, contacts are polarised to prevent incorrect connection.
Thermocouple Type	R/S IEC
Connector Type	Socket
Connector Size	Standard
Colour	Orange
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Standard Connector