

JIS Tipo K 0.7kg Botn de traccin Magnet Termopar, PFA aislado con tapn miniatura



[Click here to view product](#)

Product Code	XE-0550
General Description	Ideal for ferrous metal surface measurement
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	1/2" diameter button type, 0.7kg Pull
Cable Length	0.5m, 1m, 2m
Cable Type	PFA insulated twin-twist wires
Cable Termination	Miniature Fitted Plug
Max. Temperature	+250°C (plug 200°C)
Standards Met	JIS (Japanese Industrial Standards)

Colour	Positive: Red / Negative: White
--------	---------------------------------