


Termopar de banda magntica tipo K ANSI, aislado con PFA con enchufe en miniatura - 1M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

Product Code	XE-0579-001
General Description	Ideal for ferrous metal surface measurement
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Magnet	50 x 25 x 1mm magnetic strip
Cable Length	1M
Cable Type	PFA insulated wires, twin twisted
Termination Type	Miniature Plug
Standards Met	ANSI MC96.1
Colour	Positive: Red / Negative: Yellow
Series	Magnet Thermocouple