

**Termopar de pernos ANSI, acero inoxidable de
fibra de vidrio sobre cable trenzado - Tipo K -
Rosca M12 - Perno de 13 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

Product Code	XE-0588-001
Thermocouple Type	K
Thermocouple Junction	Bolts have grounded junctions
Cable Length	2M
Cable Type	Glassfibre stainless steel over braided cable
Cable Termination	Bare tails
Max. Temperature	250°C
Standards Met	ANSI MC96.1
Bolt thread/pitch	M12 x 1mm
Colour	Positive: Red / Negative: Yellow
Length	13mm Bolt
Series	Bolt Thermocouple