

Termopar de imn de botn de traccin tipo K de 1,9 kg ANSI, cable aislado de PFA con terminacin de trenza de acero inoxidable en enchufe en miniatura - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-0621-001
General Description	1.9kg Pull Button Magnet Thermocouple
Thermocouple Type	K
Magnet	19.1mm Diameter / 12.7mm height
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel overbraid 7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	ANSI
Series	Magnet Thermocouple