

JIS Tipo K 1.9kg Botn de tirn Magnet Termopar, Cable aislado PFA con Acero Inoxidable Sobre-Trenza Terminando en Miniatura o Enchufe Estndar



[Click here to view product](#)

Product Code	XE-0624
General Description	1.9kg Pull Button Magnet Thermocouple
Thermocouple Type	K
Magnet	19.1mm Diameter / 12.7mm height - Max pull: 1.9 KG
Cable Length	1.0m or 2.0m
Cable Type	PFA Insulated Cable with stainless steel overbraid 7/0.2mm
Termination Type	Miniature or Standard Pulg
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	JIS