

**Termopar conectado a tierra fabricado en tubo de  
acero inoxidable, acero inoxidable de fibra de  
vidrio Cable trenzado - Tipo K JIS - 4,76 mm de  
dimetro x 100 mm de longitud**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Características técnicas

**Specifications**

<b>Product Code</b>	XE-0653-001
<b>Sensor Type</b>	316 stainless steel tube in various sizes
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded at tip for fast response
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre stainless steel overbraided cable
<b>Cable Termination</b>	Bare Wire Tails
<b>Max. Temperature</b>	350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Dimensions</b>	4.76mm Diameter x 100mm Length