

Termopar de superficie de tubera ANSI, trenza de acero inoxidable de fibra de vidrio - Tipo K para adaptarse a un diametro de tubera de 30-45 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-0656-001
General Description	Stainless steel pipe clamp sensor with adjustable diameter
Thermocouple Type	K
Thermocouple Junction	Grounded tip for fast response
Cable Length	2M
Cable Type	Glassfibre stainless steel braided extension cable
Core / Strands	7/0.2mm
Termination Type	Bare Tails
Max. Temperature	+350°C
Standards Met	ANSI MC96.1
To Fit Pipe Diameters	30-45mm
Series	Pipe Surface Thermocouple