



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar de imn de traccin (herradura) tipo K ANSI tipo K, cable aislado con PFA con terminacin de trenzado de acero inoxidable en enchufe en miniatura - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-0667-001
General Description	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
Cable Length	2M
Cable Type	PFA Teflon® insulated with stainless steel over-braid
Core / Strands	7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1
Series	Magnet Thermocouple