

**Termopar magnético de potencia de tracción
(herradura) tipo K ANSI, cable aislado con PFA
con terminación de trenza de acero inoxidable en
enchufe en miniatura - 3M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-0663-001
General Description	Rapid, reliable temperature sensing on ferrous surfaces
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
Cable Length	3M
Cable Type	PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1
Series	Magnet Thermocouple