



## Conector de termopar rectangular estndar Socket de fascia IS-K-FF Tipo K IEC



**Standard Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar rectangular estndar Socket de fascia IS-K-FF Tipo K IEC**

Pinza de nylon para fijar rpido y sencillo al panel. Todos los zcalos de montaje de la fascia estn polarizados mediante pines de tamao desigual

**Especificaciones****Specifications**

<b>Product Code</b>	XE-1407-001
<b>General Description</b>	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
<b>Thermocouple Type</b>	K IEC
<b>Connector Type</b>	Fascia Mounting Socket
<b>Connector Size</b>	Standard
<b>Colour</b>	Green
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS
<b>Series</b>	Standard Connector