

**Termopar de hoja ANSI, cable aislado con fibra de vidrio con trenzado de acero inoxidable que termina en enchufe en miniatura - 1 m de longitud**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Características técnicas****Specifications**

|                              |                                                          |
|------------------------------|----------------------------------------------------------|
| <b>Product Code</b>          | XE-0686-001                                              |
| <b>General Description</b>   | Thermocouple mounted on a stainless steel shim           |
| <b>Thermocouple Type</b>     | K                                                        |
| <b>Thermocouple Junction</b> | Junction grounded at tip for fast response               |
| <b>Cable Length</b>          | 1M                                                       |
| <b>Cable Type</b>            | Glassfibre insulated lead with stainless steel overbraid |
| <b>Core / Strands</b>        | 7/0.2mm                                                  |
| <b>Termination Type</b>      | Miniature Plug                                           |
| <b>Max. Temperature</b>      | 350°C                                                    |
| <b>Min. Temperature</b>      | -50°C                                                    |
| <b>Standards Met</b>         | ANSI MC96.1                                              |
| <b>Dimensions</b>            | 25x13mm shim (0.2mm thick)                               |
| <b>Series</b>                | Leaf Thermocouple                                        |