

Termopar de hoja ANSI, cable aislado con fibra de vidrio con trenzado de acero inoxidable que termina en enchufe en miniatura - 2 m de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-0687-001
General Description	Thermocouple mounted on a stainless steel shim
Thermocouple Type	K
Thermocouple Junction	Junction grounded at tip for fast response
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	350°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1
Dimensions	25x13mm shim (0.2mm thick)
Series	Leaf Thermocouple