

Sonda de termopar de uso general ANSI, cable aislado con fibra de vidrio con trenza de acero inoxidable - Terminación tipo K en enchufe en miniatura - longitud de 1 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockists of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-0762-001
General Description	Rugged thermocouple with stainless steel sheath
Thermocouple Type	K
Thermocouple Junction	Grounded hot junction for fast response
Probe Diameter	4.5mm
Probe Length	150mm
Cable Length	1M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	ANSI MC96.1
Series	Fabricated Tube Thermocouple