



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar de boquilla ANSI, cable aislado con fibra de vidrio con trenza de acero inoxidable - Tipo K que termina en colas desnudas - 2 m de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-0870-001
General Description	Thermocouple housed in a stainless steel threaded bolt
Thermocouple Type	K
Thermocouple Junction	Grounded hot junction for fast response, M6x1mm thread
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Bare tails
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	ANSI MC96.1
Series	Nozzle Thermocouple