



Conector rectangular estndar de termopar y toma de fascia FSTC-CU-FF Tipo Cu



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

<p><resistente>Conector rectangular estndar de termopar Socket de fascia FSTC-CU-FF Tipo CU</fuerte></p><p>Clip de nailon para fijar rpido y sencillo al panel. Todos los zcalos de montaje de la fascia estn polarizados mediante pines de tamao desigual</p><p>><fuerte>Especificaciones</fuerte></envergadura></p>

Specifications

Product Code	XE-1677-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	White
Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS