




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **Sonda de termopar de uso general IEC, cable aislado con fibra de vidrio con trenzado de acero inoxidable - Tipo K - Terminación de enchufe en miniatura - 2M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Características técnicas

**Specifications**

<b>Product Code</b>	XE-0755-001
<b>General Description</b>	Rugged thermocouple with stainless steel sheath
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded hot junction for fast response
<b>Probe Diameter</b>	4.5mm
<b>Probe Length</b>	150mm
<b>Cable Length</b>	2M
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Miniature Plug (with cable clamps)
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	IEC