

**Sonda de termopar de uso general IEC, cable  
aislado con fibra de vidrio con trenzado de acero  
inoxidable - Tipo J - Terminacin de colas  
desnudas - 1M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Características técnicas****Specifications**

|                              |                                                          |
|------------------------------|----------------------------------------------------------|
| <b>Product Code</b>          | XE-0761-001                                              |
| <b>General Description</b>   | Rugged thermocouple with stainless steel sheath          |
| <b>Thermocouple Type</b>     | J                                                        |
| <b>Thermocouple Junction</b> | Grounded hot junction for fast response                  |
| <b>Probe Diameter</b>        | 4.5mm                                                    |
| <b>Probe Length</b>          | 150mm                                                    |
| <b>Cable Length</b>          | 1M                                                       |
| <b>Cable Type</b>            | Glassfibre insulated lead with stainless steel overbraid |
| <b>Core / Strands</b>        | 7/0.2mm                                                  |
| <b>Termination Type</b>      | Bare tails                                               |
| <b>Max. Temperature</b>      | +350°C                                                   |
| <b>Min. Temperature</b>      | -60°C                                                    |
| <b>Standards Met</b>         | IEC                                                      |