

Sonda de termopar de uso general IEC, cable aislado de fibra de vidrio con trenzado de acero inoxidable - Tipo J - Terminacin de enchufe estndar - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas

Specifications

Product Code	XE-5628-001
General Description	Rugged thermocouple with stainless steel sheath
Thermocouple Type	J
Thermocouple Junction	Grounded hot junction for fast response
Probe Diameter	4.5mm
Probe Length	150mm
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Standard Plug (with cable clamps)
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	IEC