

**Termopar de boquilla IEC, cable aislado de fibra
de vidrio con trenza de acero inoxidable - Tipo K -
Terminacin de enchufe estndar - 2M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

Product Code	XE-5626-001
General Description	Thermocouple housed in a stainless steel threaded bolt
Thermocouple Type	K
Thermocouple Junction	Grounded hot junction for fast response, M6x1mm thread
Probe Diameter	4.5mm
Probe Length	150mm
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Cable Termination	Standard Plug
Termination Type	Grounded hot junction for fast response, M6x1mm thread
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	IEC
Connector Type	Plug
Connector Size	Standard
Dimensions	25x13mm shim (0.2mm thick)
Thread	M6 x 1mm thread