



## Aceitunas PTFE



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Aceitunas PTFE. Disponible para adaptarse a diámetros de sonda que van desde 1,5 mm a 8,0 mm

## Specifications

<b>General Description</b>	Olives
<b>Material</b>	PTFE
<b>Probe Diameter</b>	Suitable for probe diameters between 1.5mm & 8.0mm