


**Termopar de unin expuesta IEC de doble torsin  
con aislamiento de PFA con colas desnudas -  
Tipo J - 1/0,5 mm - Longitud de cable de 1 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Termopar tipo K, J o T. Respuesta rpida, unin expuesta soldada disponible en la especificacin de cable de Clase 1 de acuerdo con BS EN 60584-3:2008 con terminacin de cola desnuda.

Ideal para aplicaciones de prueba y desarrollo

### Caractersticas tcnicas

## Specifications

|                     |                                         |
|---------------------|-----------------------------------------|
| Product Code        | XE-7135-001                             |
| General Description | Fast response, welded exposed junction  |
| Sensor Type         | Exposed Junction                        |
| Thermocouple Type   | J IEC                                   |
| Cable Length        | 1M                                      |
| Cable Type          | Teflon® PFA insulated twin-twisted lead |
| Core / Strands      | 1/0.5mm                                 |
| Cable Termination   | Bare tails                              |
| Max. Temperature    | 250°C                                   |
| Min. Temperature    | -75°C                                   |
| Tolerance           | Class 1                                 |
| Standards Met       | IEC-584-3                               |
| Series              | Exposed Junctions                       |