

Termopar de unin expuesta IEC de doble trenzado con aislamiento de PTFE con colas desnudas - Tipo T - Longitud de cable de 1 m

1000 Ohm 1000 Ohm 1000 Ohm 1000 Ohm 1000 Ohm
1000 Ohm 1000 Ohm 1000 Ohm 1000 Ohm 1000 Ohm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Termopar tipo K, J o T. Respuesta rpida, unin expuesta soldada disponible en la especificacin de cable de Clase 1 de acuerdo con BS EN 60584-3:2008 con terminacin de cola desnuda.

Ideal para aplicaciones de prueba y desarrollo.

Caractersticas tcnicas

Specifications

Product Code	XE-3520-001
General Description	Fast response, welded exposed Junction
Sensor Type	Exposed Junction
Cable Length	1M
Core / Strands	1/0.2mm
Cable Termination	Bare tails
Max. Temperature	+260°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC-584-3