

**Sondas de temperatura de termopar de rosca IEC
- Tipo K - 1/4" BSPP - 150 mm de longitud x 3 mm
de diametro**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sensores de temperatura tipo K atornillados con un rango de medicin de -75 a +250 C y cable conductor de par plano PFA de 1 metro.

Para montaje directo en recipientes, tuberas, tanques o conductos para la medicin de temperatura de liquidos o gases.

Aplicaciones tpicas

Control de procesos industriales

Tuberas o conductos

Recipientes y tanques

Caractersticas tcnicas

Specifications

Product Code	XE-7334-001
General Description	Screw in Type K thermocouple
Thermocouple Type	K IEC
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	150mm
Sheath Diameter (mm)	3mm
Process Connection options	1/4"BSPP
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	IEC
Number of Wires	2