

## **Sondas de temperatura de termopar de rosca IEC - Tipo K - 1/4"BSPP - 200 mm de longitud x 3 mm de diametro**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sensores de temperatura tipo K atornillados con un rango de medicin de -75 a +250 C y cable conductor de par plano PFA de 1 metro.

Para montaje directo en recipientes, tuberías, tanques o conductos para la medicin de temperatura de líquidos o gases.

Aplicaciones tpicas

Control de procesos industriales

Tuberías o conductos

Recipientes y tanques

#### **Caractersticas tcnicas**

## **Specifications**

|                                   |                              |
|-----------------------------------|------------------------------|
| <b>Product Code</b>               | XE-7335-001                  |
| <b>General Description</b>        | Screw in Type K thermocouple |
| <b>Thermocouple Type</b>          | K IEC                        |
| <b>Cable Length</b>               | 1M                           |
| <b>Cable Type</b>                 | PFA                          |
| <b>Termination Type</b>           | Stripped Wires               |
| <b>Sheath Length (mm)</b>         | 200mm                        |
| <b>Sheath Diameter (mm)</b>       | 3mm                          |
| <b>Process Connection options</b> | 1/4"BSPP                     |
| <b>Max. Temperature</b>           | 250°C                        |
| <b>Min. Temperature</b>           | -75°C                        |
| <b>Standards Met</b>              | IEC                          |
| <b>Number of Wires</b>            | 2                            |