

Panel de aluminio para conector de termopar con enchufes estándar IEC tipo N - 20 vas



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Paneles de conectores estndar de termopar preensamblados

Un conjunto estndar de panel horizontal preensamblado diseado para usarse con termopares o cables con tapones estndar. Robusto y perfecto para montarse en salas de control, bancadas de pruebas, aplicaciones de laboratorio, salas de control y ms. Admiten cdigos de calibracin como K, J, T & N con enchufes de cable rpidos y pueden acomodar circuitos de 1-24 vas, con un rango mximo de temperatura de hasta +120C.

Suministrado previamente equipado con

[vasos estndar de fascia](#)

Circuito de 20 vas

Calibracin tipo N

Temperatura mxima 120C

Aplicaciones tpicas

OEMs

Constructores de paneles de instrumentacin

Salas de control

Bancos de pruebas y patines

Specifications

Product Code	XE-7567-001
Thermocouple Type	N IEC
Number of Ways	20 Way
Connector Size	Standard
Connector Type	Fascia Mounting Socket
Colour	Pink
Dimensions	H111mm x W221mm
Max. Temperature	120°C
+Positive Contact	Nickel / Chromium / Silicon



-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS