

Panel de aluminio para conector termopar con tomas miniatura ANSI tipo T - 18 vas



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Paneles de conectores de termopar miniatura preensamblados

Un conjunto de panel horizontal en miniatura preensamblado diseñado para usarse con termopares o cables con tapones en miniatura. Robusto y perfecto para montarse en salas de control, bancadas de pruebas, aplicaciones de laboratorio, salas de control y ms. Admiten cdigos de calibracin como K, J, T con conectores rpidos y pueden acomodar circuitos de 1 a 24 vas, con un rango mximo de temperatura de hasta +120C.

Suministrado previamente equipado con

[vasos de fascia miniatura](#)

Circuito de 18 vas

Calibracin ANSI tipo T

Temperatura mxima 120C

Aplicaciones tpicas

OEMs

Constructores de paneles de instrumentacin

Salas de control

Bancos de pruebas y patines

Specifications

Product Code	XE-7553-001
General Description	A selection of miniature IEC pre-assembled thermocouple panels, 1-24 way, to suit your specific application. All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	T ANSI
Number of Ways	18 Way
Connector Size	Miniature
Connector Type	Fascia Mounting Socket
Colour	Blue



Dimensions	H75mm x W160mm
Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS