

Panel de aluminio con conector de termopar con enchufes estándar ANSI tipo T - 3 vas



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Paneles de conectores estndar de termopar preensamblados

Un conjunto estndar de panel horizontal preensamblado diseado para usarse con termopares o cables con tapones estndar. Robusto y perfecto para montarse en salas de control, bancadas de pruebas, aplicaciones de laboratorio, salas de control y ms. Admiten cdigos de calibracin como K, J, T con conectores rpidos y pueden acomodar circuitos de 1 a 24 vas, con un rango mximo de temperatura de hasta +120C.

Suministrado previamente equipado con

[vasos estndar de fascia](#)

Circuito de 3 vas

Calibracin ANSI tipo T

Temperatura mxima 120C

Aplicaciones tpicas

OEMs

Constructores de paneles de instrumentacin

Salas de control

Bancos de pruebas y patines

Specifications

Product Code	XE-7519-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel.
Thermocouple Type	T ANSI
Number of Ways	3 Way
Connector Size	Standard
Connector Type	Fascia Mounting Socket
Colour	Blue
Dimensions	H66mm x W88mm



Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS