

Panel de aluminio para conector termopar con tomas miniaturas JIS Tipo K - 16 direcciones



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Paneles de conectores de termopar miniatura preensamblados

Un conjunto de panel horizontal en miniatura preensamblado diseñado para usarse con termopares o cables con tapones en miniatura. Robusto y perfecto para montarse en salas de control, bancadas de pruebas, aplicaciones de laboratorio, salas de control y ms. Admiten el código de calibración K con cables de cable rápido y pueden acomodar circuitos de 1-24 vas, con un rango máximo de temperatura de hasta +120C.

Suministrado previamente equipado con

[vasos de fascia miniatura](#)

Circuito de 16 vas

Calibración JIS tipo K

Temperatura máxima 120C

Aplicaciones típicas

OEMs

Constructores de paneles de instrumentación

Salas de control

Bancos de pruebas y patines

Specifications

Product Code	XE-7604-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	K JIS
Number of Ways	16 Way
Connector Size	Miniature
Connector Type	Fascia Mounting Socket
Colour	Blue
Dimensions	H75mm x W145mm
Max. Temperature	120°C



+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS