

Panel de aluminio para conector termopar con enchufes estndar JIS tipo K - 10 vas



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Paneles de conectores estndar de termopar preensamblados

Un conjunto estndar de panel horizontal preensamblado diseado para usarse con termopares o cables con tapones estndar. Robusto y perfecto para montarse en salas de control, bancadas de pruebas, aplicaciones de laboratorio, salas de control y ms. Admiten el cdigo de calibracin K con zcalos de cable rpido y pueden acomodar circuitos de 1-24 vas, con un rango mximo de temperatura de hasta +120C.

Suministrado previamente equipado con

[vasos estndar de fascia](#)

Circuito de 10 vas

Calibracin tipo K

Temperatura mxima 120C

Aplicaciones tpicas

OEMs

Constructores de paneles de instrumentacin

Salas de control

Bancos de pruebas y patines

Specifications

Product Code	XE-7589-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel.
Thermocouple Type	K JIS
Number of Ways	10 Way
Connector Size	Standard
Connector Type	Fascia Mounting Socket
Colour	Blue
Dimensions	H66mm x W221mm



Max. Temperature	120°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS