

Termopares de tubo fabricados por IEC con sello de engarce y cable conductor - Tipo J - Dimetro de la funda de 4,5 mm - Longitud de la funda de 50 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Aplicaciones tpicas
Control de Procesos Industriales
Tuberias o conductos
Hornos y calefaccin
Propsito general

Caractersticas tcnicas**Specifications**

Product Code	XE-7916-001
General Description	Fabricated IEC Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead
Thermocouple Type	J
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	50mm
Sheath Diameter (mm)	4.5mm
Max. Temperature	250°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC
Number of Wires	2
Series	Fabricated Tube Thermocouple