

Termopares de tubo fabricados por ANSI con sello de engarzado y alambre conductor - Tipo K - 4,5 mm x 250 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Aplicaciones tpicas

Control de procesos industriales

Tuberas o conductos

Hornos y calefaccin

Uso general

Caractersticas tcnicas**Specifications**

Product Code	XE-7943-001
General Description	Fabricated Type K Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead
Thermocouple Type	K ANSI
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	50, 100, 150, 250
Sheath Diameter (mm)	3.0, 4.5, 6.0
Max. Temperature	250°C
Min. Temperature	75°C
Standards Met	Class 1
Number of Wires	2
Series	Fabricated Tube Thermocouple