


## Conector de Termopar Miniatura Tipo Toma en Lnea AM-R/S-FS R/S ANSI



**Miniature in-line Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar miniatura Enchufe en Inea AM-R/S-FS Tipo R/S ANSI**

Ideal para uso general, todos los contactos estn polarizados para asegurar una conexin correcta. Estos conectores aceptan cables de termopar de hasta 4 mm de dimetro.

**Especificaciones****Specifications**

<b>Product Code</b>	XE-1312-001
<b>General Description</b>	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
<b>Thermocouple Type</b>	R/S ANSI
<b>Connector Type</b>	Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Green
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS
<b>Series</b>	Miniature Connector