



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar de unin expuesta IEC con doble trenzado aislado PTFE y colas desnudas - Tipo K - 300 mm de longitud de cable



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Termopar tipo K, J o T. Respuesta rpida, unin soldada expuesta, disponible en la especificacin de cable Clase 1 conforme a la BS EN 60584-3:2008 con terminacin de cola desnuda.

Ideal para aplicaciones de prueba y desarrollo.

Especificaciones

Specifications

Product Code	XE-0424-001
General Description	Fast response, welded exposed Junction
Sensor Type	Exposed Junction
Cable Length	300mm
Core / Strands	1/0.2mm
Cable Termination	Bare tails
Max. Temperature	+260°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC-584-3
Series	Exposed Junctions